

Validate XP Software




Expedite Control Software For Windows XP

Validate XP Software offers the full range of control and monitoring of synthesis operations on the Expedite synthesizer while taking advantage of all the modern features of Windows XP.



Validate XP offers flexible communication options to the Expedite through ISA or PCI expansion bus or standard 9 pin RS232



Validate XP

- Controls** *Expedite, Expedite-MOSS, and Expedite-ITEM configurations.*
- Displays** *active synthesizer status and saves all synthesis data to a Microsoft Access database format.*
- Provides** *several built-in report formats including the synthesis data report as well as a performance and trityl reports.*
- Utilizes** *the standard Expedite protocol files and protocol editor so that you can continue to use the same proven process developed for your synthesizer.*

Genomic Technologies, Inc.

12 Walnut Hill Road
Millis, MA 02054
617.901.6178
www.GenomicTechnologies.com
Info07@GenomicTechnologies.com

Providing Cost Effective Solutions for the High Throughput Service Laboratory